



INNOVATIVE SLIDE FORMATTING



Samples of 12, 48 & 192 arrays

Design	Volume per well
12 array	40ul
16 array	60ul
48 array	6-7ul
60 array	5ul
192 array	1-1.5ul

- imprint patterns onto slides (any design/pattern can be imprinted on slides)
- high throughput sample processing
- parallel loading of samples



REPLACEMENT STAMPS

- 12 Arrays
CAT# WSP12-R
- 16 Arrays
CAT# WSP16-R
- 48 Arrays
CAT# WSP48-R
- 60 Arrays
CAT# WSP60-R
- 192 Arrays
CAT# WSP192-R

Slidelmprinter

12 Arrays
CAT# WSP12-1

16 Arrays
CAT# WSP16-1

48 Arrays
CAT# WSP48-1

60 Arrays
CAT# WSP60-1

192 Arrays
CAT# WSP192-1

Custom Design
CAT# CUST-1

Hydrophobic Wax - Low Temperature
CAT# HWAX

DNA Wax - High Temperature
CAT# DWAX

Test Slides - 72/pkg
CAT# TSLD



Slide Imprinter



665 Third Street, Suite 240 | San Francisco, CA 94107 | USA
415-247-8760 (tel) | 415-247-8765 (fax)
www.gelcompany.com




What is Slidelmprinter


The Slidelmprinter is a much more effective and productive way to partition slides used in laboratory and microarray analysis. This instrument is ideal for high throughput screening where multiple assays need to be run on a single examination slide with minimal risk of contamination between assays.

Traditionally, rubber gaskets (with or without adhesive), hydrophobic markers, or manufactured (teflon) barriers have been used to segregate slide areas. Gaskets without adhesive may not seal well. Gaskets with adhesive may also leak as the adhesive can be soluble in some solvents and interfere with the experiments. Gaskets also need to be removed before scanning which can cause the entire slide coating to peel. Hydrophobic markers (PAP pens) are very tedious and messy to apply, can be soluble in certain solvents, and the resulting slides are often irreproducible and inaccurate. Lastly, manufactured barriers are expensive and not flexible to design.

The Slidelmprinter provides a convenient way to define multiple “wells” on a single microscope slide to create separated region for sampling. This is done through imprinting a thin, inert, hydrophobic



- 1) Slide inserted into slide holder
- 2) Stamp with the specific design
- 3) Handle to lift stamp up and down




Sample loading onto a 48 well imprinted slide.

16 array test slides


How does Slidelmprinter work

film line onto the surface of a slide through the use of a designed “stamp.” The stamp is immersed in a bath. A slide is inserted into the slide holder above the bath. The stamp is lifted to meet the slide, imprinting barrier lines. Since the deposits are thin, the resulting slide with lines can be inserted directly into a slide scanner without any modifications. The coating, being highly hydrophobic, does not dissolve in most solvent systems, and will stay on the slide through multiple wash steps.

The Slidelmprinter offers flexibility in slide partition design, adhesion with any slide surface chemistry, stability in most solvents through repeated washes, is chemically inert to testing, and the barrier does not need to be removed prior to scanning.



Insert slide into slide holder prior to imprinting. (active side faced down)



Five stamp configurations available:
12, 16, 48, 60 & 192 arrays

The Slidelmprinter offers flexibility in slide partition design, adhesion with any slide surface chemistry, stability in most solvents through repeated washes, is chemically inert to testing, and the barrier does not need to be removed prior to scanning.