




# Omega Fluor™ and Omega Fluor Plus™

Gel Documentation Systems

User Manual



## Table of Contents

Section	Page
<b>Omega Fluor</b>	
Introduction	3
Components and Parts	3
Omega Fluor System Specifications	4
System Placement	5
Setup Illustrations	5
Software Installation	6
<b>Omega Fluor Plus</b>	
Introduction	12
Components and Parts	12
Omega Fluor Plus System Specifications	13
System Placement	14
Setup Illustrations	14
<b>Omega Fluor and Omega Fluor Plus</b>	
Taking an Image	16
<b>Omega Fluor</b>	
Taking an Image on the Omega Fluor	17
Selecting Your Light Source	17
Capturing Your Image	18
<b>Omega Fluor Plus</b>	
Taking an Image on the Omega Fluor Plus	20
Selecting Your Light Source	20
Capturing Your Image	22
<b>Omega Fluor and Omega Fluor Plus</b>	
Overview of Software Features	23
Omega Fluor UV Cutoff Switch	27
Common Imaging Questions	28
Appendix A—Ordering Information	29
Appendix B—Routine Maintenance	30
Appendix C—Regulatory	33
Appendix D—Mitsubishi Electric P-95D Printer Driver Installation	36
<b>Omega Fluor</b>	
Appendix E—System Requirements	39
<b>Omega Fluor Plus</b>	
Appendix F—Installing or Changing Emission Filters	40
Appendix G—USB Installation Guide	41
Appendix H—Basic Tablet Operations	43

## Introduction

Congratulations on purchasing your Aplegen Omega Fluor Gel Documentation System. The Omega Fluor is a powerful yet simple tool for gel documentation and generating publication quality, 16bit tiff images. The Omega Fluor comes completely assembled for quick startup. It includes Omega Fluor Acquisition Software which is compatible with Microsoft Windows XP and Vista, Windows 7 (32 and 64 bit), and Windows 8 (32bit only). It has a clean user interface and simple tools for annotations and contrast adjustments. For more in depth analysis, the data can be exported to a gel analysis software.

## Components and Parts

The Omega Fluor system includes the following components:

- **The cabinet**—The cabinet is a light tight imaging station. It includes a 302nm UV transilluminator on a pull out tray. Samples can be seen directly through the UV protected viewport. The viewport includes an orange filter for improved visual contrast. The cabinet also contains an EPI white light.
- **The camera**—The Omega Fluor camera is a 5MP and is pre-installed in the cabinet. The 590/50nm orange filter is pre-assembled on the end of the lens
- **The software**—The Omega Fluor Acquisition Software is included on the USB Drive and can be loaded on any Windows™ compatible computer.



## Omega Fluor System Specifications

Specifications	
Camera	<ul style="list-style-type: none"> <li>• 5MP for high resolution images</li> <li>• USB connection for simple set up</li> <li>• Fast refresh rate for live imaging</li> <li>• High quality images ideal for downstream analysis</li> </ul>
Lens	8mm F1.4 Lens
Field of View	20 cm by 24 cm
Cabinet	<ul style="list-style-type: none"> <li>• 302 nm pull out UV transilluminator (or optional 365 nm)</li> <li>• UV viewing port</li> <li>• EPI white lights</li> </ul>
Emission Filter	590 nm
Applications	<ul style="list-style-type: none"> <li>• Fluorescent gel imaging</li> <li>• Visible gel imaging</li> <li>• Gel documentation</li> </ul>
Certifications	CE; CSA; UL compliant
Product Footprint	34.6cm x 31.1cm x 68.6cm
Computer Requirements	<p>Desktop Computers</p> <ul style="list-style-type: none"> <li>• Supported OS: Windows 8 32 bit; Windows 7/Vista 32 or 64 bit, Windows XP (with Microsoft .NET Framework 3.5)</li> <li>• Minimum hardware requirements: 1.4 GHz processor speed, 2GB RAM, 16 GB Free Hard Disk Space, 2 USB (camera and printer)</li> </ul> <p>Laptop Computers</p> <ul style="list-style-type: none"> <li>• Supported OS: Windows 8 32 bit; Windows 7/Vista 32 or 64 bit, Windows XP (with Microsoft .NET Framework 3.5)</li> <li>• Minimum hardware requirements: 1.4 GHz processor speed, 2GB RAM, 16 GB Free Hard Disk Space, 2 USB (camera and printer)</li> </ul>

## System Placement

As with all electrical instruments, the Omega Fluor imaging system should be located away from water, solvents, or corrosive materials, on a flat and stable surface with adequate clearance on all sides. The top of the system should have at least 10 cm clearance to allow sufficient air flow around the camera head.

The system is intended for indoor use with the following ambient conditions:

- Altitude up to 2000 m;
- Temperature 5 °C to 40 °C;
- Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C;

Further, the system should be placed away from interfering electrical signals and magnetic fields. If possible, a dedicated electrical outlet should be used to eliminate electrical interference from other instrumentation in your laboratory.

## Setup Illustrations

- Begin by placing the Omega Fluor on a flat, level surface with adequate clearance on all sides.
- Plug the USB camera connection into the computer.
- Connect the power supply and the AC adapter cord.
- Plug the cord into a power outlet.
- **Aplegen recommends a surge protecting powerstrip to protect against potential damage from power surges.**
- Install the Omega Fluor Software onto the computer which will run the instrument.



## Software Installation




The Omega Fluor Acquisition Software is provided on the included USB drive. Installation should be completed by a user with administrator rights. A complete list of system requirements can be found in Appendix D.

### Windows 8, 32 Bit Installation


To install the software on a Windows 8, 32 bit operation system you will need to disable the Secure Boot feature prior to installing the software and drivers. The Secure Boot feature is an advanced security feature new to the Windows 8 operating system.

To disable the Secure Boot:


- Open the Charms Bar, swipe with your finger from the far right side of the screen, to the left.
- Select **Settings**.
- Select **Change PC Settings**.
- Select **General**.
- Scroll down the list, Under Advanced Startup select **Restart Now**. The system will restart.




- Select **Troubleshoot**.




- Select **Advanced Options**.




- Select **Startup Settings**.



- Select **Restart**.




- Use the USB keyboard to select number 7 (hit **F7** or the number 7 key). The system will restart. You are now ready to install the printer driver.
- To start, plug in the printer USB.
- Turn on the printer and follow the directions that appear on the screen.







With the Secure Boot disabled, plug in the USB.




- Select the **Omega Fluor Setup Program** from the launch window.



- Follow the directions on the screen.



- You are now ready to use your Omega Fluor.






## Windows 7 Software Installation

- To complete Windows 7 software installation, plug in the USB drive.



- Follow the instructions on the screen.



A pop up may appear that says the Omega Fluor Camera has not passed Windows Logo testing. The camera driver is a standard camera driver and has been extensively tested to ensure compatibility. In addition, it is currently in process of Windows Logo verification. Click **Continue Anyway** to complete installation.




## Windows XP Software Installation


- To install the software on a Windows XP machine, plug in the USB.



- When prompted, select **Omega Fluor Setup Program** from the launch window.



- Follow the instructions on the screen.







You may also navigate directly to the OmegaFluorGDS-Setup.exe through the Computer menu.

- Select **My Computer** from the Application Menu
- Select the **Omega Fluor (G:)** drive
- Select the **OmegaFluorGDS-Setup.exe** and double click to run.

The first time the camera is connected the Found New Hardware Wizard will pop up to install the camera driver.

- To complete the installation select **No not this time**.
- Select **Install Software Automatically (Recommended)**.
- When prompted click **Finish** to complete installation of the camera driver.

A pop up may appear that says the Omega Fluor Camera has not passed Windows Logo testing. The camera driver is a standard camera driver and has been extensively tested to ensure compatibility. In addition, it is currently in process of Windows Logo verification. Click **Continue Anyway** to complete installation.



## Introduction

Congratulations on purchasing your Aplegen Omega Fluor Plus gel documentation system. The Omega Fluor Plus is a powerful yet simple tool for gel documentation and generating publication quality, 16bit tiff images. The Omega Fluor Plus comes completely assembled for quick startup, and includes Omega Fluor Acquisition Software. It has a clean user interface and simple tools for annotations and contrast adjustments. For more in depth analysis, the data can be exported to a gel analysis software.

## Components and Parts

The Omega Fluor Plus system includes the following components:


- **The cabinet**—The cabinet is a light tight imaging station. It includes a 302nm UV transilluminator on a pull out tray. The cabinet also contains an EPI white and EPI Blue light.



- **The camera**—the Omega Fluor Plus camera is a 5MP, and comes installed on the cabinet. There is a 3 position filter wheel with an Orange 590/50 nm filter standard.



- **The software**—The Omega Fluor Plus Acquisition Software is included on the USB Drive and can be loaded on any Windows™ compatible computer.



## Omega Fluor Plus System Specifications

Specifications	
Camera	<ul style="list-style-type: none"><li>• 5MP for high resolution images</li><li>• USB connection for simple set up</li><li>• Fast refresh rate for live imaging</li><li>• High quality images ideal for downstream analysis</li></ul>
Lens	8mm F1.4 Lens
Field of View	20 cm by 24 cm
Cabinet	<ul style="list-style-type: none"><li>• 302 nm pull out UV transilluminator (or optional 365 nm)</li><li>• EPI white lights</li></ul>
Emission Filter	590 nm
Applications	<ul style="list-style-type: none"><li>• Fluorescent gel imaging</li><li>• Visible gel imaging</li><li>• Gel documentation</li></ul>
Certifications	CE; CSA; UL compliant
Product Footprint	34.6cm x 31.1cm x 68.6cm

## System Placement

As with all electrical instruments, the Omega Fluor Plus imaging system should be located away from water, solvents, or corrosive materials, on a flat and stable surface with adequate clearance on all sides. The top of the system should have at least 10cm clearance to allow sufficient air flow around the camera head.

The system is intended for indoor use with the following ambient conditions:

- Altitude up to 2000 m;
- Temperature 5 °C to 40 °C;
- Maximum relative humidity 80 % for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C;

Further, the system should be placed away from interfering electrical signals and magnetic fields. If possible, a dedicated electrical outlet should be used to eliminate electrical interference from other instrumentation in your laboratory.

## Setup Illustrations

- Begin by placing the Omega Fluor Plus on a flat, level surface with adequate clearance on all sides.

Make sure to check inside the system and the tablet for additional packing materials and remove them prior to use.

- Connect the camera cable. You will find it just below the tablet on the right side.





- Connect the power supplies.

**Aplegen recommends a surge protecting powerstrip to protect against potential damage from power surges.**



- Turn on the power to the tablet computer. Select the Omega Fluor tile from the startup screen. You are now ready to start using your Omega Fluor Plus.



## Taking an Image

The Omega Fluor Acquisition Software allows the user to take images with their Omega Fluor system, modify the visual aspects such as contrast and saturation, make general annotations and has simple arithmetic options. The user can also open and modify previously acquired images.

### Basic Imaging involves:

1. Centering your sample
2. Selecting your lighting source
3. Capturing your image

### Basic Applications

The Omega Fluor is compatible with multiple sample types. Use the table below as a guide.

	UV Light	White Light Conversion Screen	Blue Light
Ethidium bromide	✓		
SYBR® Green	✓		✓
Coomassie blue	✓	✓	
Protein gels	✓	✓	
Visible light gels	✓	✓	
Other UV gels	✓		✓



## Taking an Image on the Omega Fluor

### Centering Your Image

Place your gel either directly on the UV transilluminator or in one of the optional sample trays available from Aplegen.

- Gel Handling Tray for UV Applications
- Gel Handling Tray for White Light Applications
- Blue Light Conversion Screen



### Selecting Your Light Source

Use this table to determine the lighting source for your application.


Dye	Excitation	Emission
Ethidium Bromide	302/365 nm	590 nm
Coomassie Blue	White Light	
Silver Stain	White Light	
SYBR® Gold	Blue Light	540 nm
SYBR® Green	Blue Light	520 nm
SYBR® Safe	Blue Light	530 nm
GelStar®	302 nm	520 nm
SYPRO® Ruby	302 nm or Blue Light	610 nm
Deep Purple™	UV or Blue Light	610 nm

## Capturing Your Image


The Omega Fluor UV transilluminator emits at a wavelength of 302nm (365nm optional). It can be controlled through the Image Capture Window. Once you have selected the appropriate light source and sample tray or conversion screen, launch the Omega Fluor Acquisition Software.




- Click **Acquire**.



- Select **UV light**.



- Click **Auto Exposure** or use the manual exposure icons to adjust the image preview.



- When using Auto Exposure the rotating marquee will appear in the bottom left of the Image Capture window.




- When the calculation is complete the marquee will change from red to green.



- If satisfied with the preview, select **Image Capture**. Your image will be captured and open in a new window.

Once the image is captured it can be further adjusted for contrast, annotations added, saved, and printed.



## Taking an Image on the Omega Fluor Plus

### Centering Your Image

Place your gel either directly on the UV transilluminator or in one of the optional sample trays available from Aplegen.

- Gel Handling Tray—UV
- Gel Handling Tray—White Light Conversion Screen



### Selecting Your Light Source

Use this table to determine the lighting source for your application.

Dye	Excitation	Emission
Ethidium Bromide	302/365 nm	590 nm
Coomassie Blue	White Light	
Silver Stain	White Light	
SYBR® Gold	Blue Light	540 nm
SYBR® Green	Blue Light	520 nm
SYBR® Safe	Blue Light	530 nm
GelStar®	302 nm	520 nm
SYPRO® Ruby	302 nm or Blue Light	610 nm
Deep Purple™	UV or Blue Light	610 nm

The Omega Fluor Plus has three lighting sources: the UV transilluminator, the EPI Blue lights, and the EPI White lights.


- Begin the acquisition by selecting your light source.




- Select your filter.



- Click the **Acquire** icon.






- The light source is controlled by the software. To turn on the light source, check the **Illumination** button.



## Capturing Your Image

- Click on the **Auto Exposure** icon for a preview of the image or use the manual exposure icons to adjust the image preview.
- When using Auto Exposure the rotating marquee will appear in the bottom left of the Image Capture window.
- When the calculation is complete the marquee will turn from red to green.
- If satisfied with the preview, select the **Image Capture** button. Your image will be captured and open in a new window.

Once the image is captured it can be further adjusted for contrast, annotations added, saved, and printed.



## Overview of Software Features

The Omega Fluor Software has a simple user interface. The top left Omega Fluor icon accesses the Application Menu, while all image acquisition and annotation functions can be accessed from the home tab. Software functions include:

### File Tab



- **Open**—Open previously saved files.



- **Save**—Save image file, or changes to files. Default file type is .tiff.



- **Save as**—Save a copy of an image file, as either a .tiff or .jpeg, or to an alternative file location.



- **Close**—Close the image file in the current tab.



- **Close all**—Closes all open image files. Unsaved images will not close without prompting the user to save or discard.

### Print



- **Print**—Print the currently open image tab.



- **Page Setup**—Change page layout and print options.



- **Print Preview**—Previews the page layout and how the image will appear on the page.



- **Create Flat**—Create flat field correction file for image modification.



- **Apply Flat**—Applies flat field correction to open image.




- **E-mail Settings**—Configure the system settings for easy export of images. Contact your network administrator if you are unsure of the SMTP server or outgoing port. The e-mail address used will appear as the sender.



- **User Manual**—Opens a PDF of the user manual.



- **About**—Provides the user with the current software version number.



### Home Tab

#### Image Menu



- **Acquire**—Click the acquire button to launch the Image Capture Window. A preview of your sample will appear.



- **Open**—Opens an image from a saved file location.



- **Save**—Opens the Save or Save As dialogue boxes.



- **Print**—Opens the Print, Print Setup, or Print Preview Dialogue boxes.



- **Crop**—Opens the Crop Image tool. To crop, select the icon. A + sign will appear in place of the mouse arrow when hovering over the image. Click and drag the mouse across the desired area. A new image will open in a new tab.



- **Resize**—An image can be resized up to 3 times its original size, or to a max size of 7776 – 5832 pixels.



- **Resize to size**—Manually enter the size in pixels of the image desired.
  - Width—Image size in the horizontal plane.
  - Height—Image size in the vertical plane.
  - Keep aspect ratio—Maintains the proportions of the original image



- **Rotate**

- Rotate Right 90°
- Rotate Left 90°
- Rotate 180°
- Rotate Arbitrary—Adjust the image angle via slide bar or by typing the desired angle in the text box.



- **Flip Horizontal**—Flip the image along the horizontal axis.

- **Flip Vertical**—Flip the image along the vertical axis.

- **Image Info**

- Displayed Information—Image Width, Image Height, Bit Depth, Exposure Time (ms), Date Taken, Taken By, Comment.
- Edit Image Information—Allows Edits to the Taken By and Comments Section of the Image Info.



- **Undo**—Undoes the last action completed by the tools in this menu.

- **Contrast Adjustments Menu**



- **Auto Contrast**—Calculates a balanced contrast ratio.



- **Negative**—Gives the inverse of the image.



- **Display Saturation**—Oversaturated pixels appear in red, under-saturated pixels appear in green.



- **Manual Contrast**—Adjust black white and gamma levels via slider or entering number into text box.

- **Zoom Menu**



- **Zoom In**—Zooms in on the image.



- **Zoom Out**—Zooms out on the image.



- **Zoom Best Fit**—Adjusts zoom to fill window. Displays % in the lower left corner of the window.



- **Panning**—Click and hold to move image right or left.



- **Annotations Menu**



- **Pointer**—This is the object selector. It allows selection of text boxes, objects, lines, etc in order to edit or delete them. For example, to change the color of a box, select the pointer, double click on the box to be changed, then select the color select tool and click on the color desired.



- **Rectangle**—Draw a box on the image.



- **Ellipse**—Draw a circle on the image.



- **Line**—Draw a line on the image.



- **Pencil**—Draw freehand.



- **Text**—Click on the text box to get a cursor. Click on the image where you want the text to appear and a text box will open. Select font and make edits to the text within the box.



- **Rotate Clockwise**—Rotates text box or other selected object clockwise.



- **Rotate Counterclockwise**—Rotates text box or other selected object counterclockwise.



- **Pen Color**—Select the color for objects or lines.



- **Pen Width**—Select the width of the line for your objects or lines.



- **Delete**—Deletes the selected object. To select an object, use the pointer, click on the desired object, then click on the Delete icon.



- **Undo**—Undoes the last action.



- **Redo**—Redoes the last action.



- **Send File Menu**—Use to send files. Email must be setup and file saved before it can be sent. See page 17 for instructions on setting up the send e-mail option.



- **The Image Capture Window**—Clicking the Acquire icon on the Home Tab launches the Image Capture Window.



- **UV Light**—Turns on the UV light for excitation of fluorescent dyes.



- **Manual Exposure**—Select an exposure time by sliding the tab along the manual exposure bar or by entering a number into the box. The range of exposure times is: 1-2394ms.



- **Auto Exposure**—Automatically determine optimal exposure time. An animated marquee appears during the calculation process. When the animation is red, the auto exposure is calculating. When the animation turns green and says “Auto exposure completed!” the system is ready to image.



- **Show Grid**—A 4x4 grid will be overlaid on the image.



- **Zoom In**—Zoom in on the desired area, the capture size is indicated on the lower left frame.



- **Zoom Out**—Zoom back out up to the original size.



- **Capture Image**—Capture the image at the selected parameters.

## Omega Fluor UV Cutoff Switch

The Omega Fluor has a UV cutoff switch to protect against accidental UV exposure when opening the cabinet door.

The cutoff switch can be overridden to allow the user to visualize or excise bands.

Before you override the UV cutoff switch make sure you are wearing the proper UV protective equipment. UV light can damage living tissue.

To override the cutoff switch:

1. Open the door and pull out the UV transilluminator door.
2. Place your sample on the UV surface.
3. Slide the UV transilluminator in and close the door.
4. Launch the software, click **Acquire** and select **UV light**. The light will automatically shut off after 5 minutes. If you need more time, select **UV light** again.



## Common Imaging Questions

### What is the difference between saving images as a .tiff or a .jpeg?

The main difference between a .tiff and a .jpeg is the amount of data processing and data compression. While the .tiff format retains information on saturation and contrast it does not compress the data. This means the .tiff format retains all of the detail in the raw image data file. The .tiff format can be saved as 8 bit or 16 bit depending on the amount of data that is collected by the sensor. The amount of data can be considerable. The .jpeg file by contrast, compresses the image into an 8 bit format. This format has reduced detail but is quickly read by image processing software, and the smaller image size allows for increased storage capacity and ease of sharing. Compression is a bit of a misnomer, as once the data has been compressed, it cannot be returned to its original state.

### I saved my image as a .jpeg and now I cannot make changes to my annotations, why?

When you save your image as a .jpeg the image both the image data and the annotation data are compressed into a single 24 bit RGB file. The annotations have functionally become part of the image. If you save the image as a .tiff the annotations are saved as a separate part of the file and can be edited by clicking on the box to select it, and making the changes as desired.

### My image seems to be fuzzy, unclear, or unfocused. How can I focus the lens?

The lens does not have a focus ring and comes pre-assembled with the correct backfocus. If the calibration is off, this can be corrected by following the procedure in Appendix B.

### I see a half moon shape on my image, bright blur, or other artifact in my image. What's wrong?

The most common cause of these artifacts are light leaks. In order to create a light tight seal around the lens there is a rubber gasket. Over time this rubber will age and cracks can form. Please follow the procedure in Appendix B, Replacing the Rubber Lens Gasket.

### When I look through the View Port I can see UV bulbs, but they don't appear in my final image. Why not?

The orange filter on the view port blocks UV and improves visual contrast. However it does not block all wavelengths of light that the bulbs emit. Therefore, the bulbs are still visible. The emission filter on the end of the Omega Fluor camera lens by contrast eliminates all emissions, except in the 590/50nm range, so the bulbs are not visible in the final image.

## Appendix A—Ordering Information

Product	Description	Part Number
Omega Fluor 302 nm configuration	Includes imaging system with cabinet, camera and trans 302 nm UV and EPI White Light sources, and Omega Fluor capture software.	
Omega Fluor 365 nm configuration	Includes imaging system with cabinet, camera and trans 365 nm UV and EPI White Light sources, and Omega Fluor capture software.	
Omega Fluor Plus 302 nm configuration	Includes imaging system with cabinet, camera and trans 302 nm UV and EPI White Light sources, and Omega Fluor Plus capture software.	
Omega Fluor Plus 365 nm configuration	Includes imaging system with cabinet, camera and trans 365 nm UV and EPI White Light sources, and Omega Fluor Plus capture software.	
Accessories		
Thermal Printer	Print Method: Thermal printing on thermal sensitive paper; Dot Density: 1280 dot/100 mm (325 dpi). USB connection.	
Thermal Paper	4 rolls	
Gel Tray	Gel Handling Tray for UV excited gels	
White Light Conversion Screen	White light conversion screen for gels requiring trans-white illumination	
Blue Light Conversion Screen	Converts UV light to blue light for use with Safe Dyes	
Desktop computer	<ul style="list-style-type: none"> <li>Operating system: Windows 8 32 bit, Windows 7 32 or 64-bit, Windows Vista 32 or 64-bit, Windows XP (with Microsoft .NET Framework 3.5).</li> <li>Minimum hardware requirements: 1.4 GHz processor speed, 1GB RAM, 16 GB Free Hard Disk Space, 2 USB (camera and printer)</li> </ul>	
Laptop Computer	<ul style="list-style-type: none"> <li>Operating system: Windows 8 32 bit, Windows 7 32 or 64-bit, Windows Vista 32 or 64 bit, Windows XP (with Microsoft .NET Framework 3.5).</li> <li>Minimum hardware requirements: 1.4 GHz processor speed, 2GB RAM, 16 GB Free Hard Disk Space, 2 USB (camera and printer)</li> </ul>	

## Appendix B—Routine Maintenance

### Cleaning the Omega Fluor or Omega Fluor Plus

Basic care is all that is required to keep your Omega Fluor or Omega Fluor Plus in great shape. Clean up any spills as they happen. Wipe down the UV glass with a soft cloth only to prevent scratching. A simple green solution can be used to clean heavier soil.

### Replacing the UV Bulbs

1. Disconnect the power to the Omega Fluor or Omega Fluor Plus.
2. Open the main door and pull out the UV transilluminator.
3. Remove the 4 screws on the two sides of the UV transilluminator.
4. Remove the top cover. Exercise caution as the top cover has a glass plate.
5. Replace the failed bulb.
6. Reinstall the UV transilluminator top cover and close the UV transilluminator drawer.

### Replacing the Rubber Lens Gasket

1. Ensure that the system is powered down.
2. Unscrew the camera attachment knob at the back of the system. Keep one hand on the camera to prevent dropping.
3. Gently remove the camera from the gasket.
4. Using a Phillips head screwdriver remove the 4 screws from the top plate, holding the gasket in place.
5. Remove and discard the rubber gasket.
6. Center the new gasket in place. There are pre-cut screw holes in the gasket, ensure they line up with the corresponding holes on the cabinet.
7. Replace the top plate and insert the 4 screws.
8. Replace the camera so that the gasket is snug and flush on the lens flashing.
9. Perform the Lens Calibration Procedure.


### Omega Fluor or Omega Fluor Plus Lens Calibration Procedure

The lens on the Omega Fluor or Omega Fluor Plus is pre-assembled on the cabinet and calibrated prior to shipment. However, it may occasional become necessary to recalibrate the lens, after changing the gasket for example. To calibrate the lens follow the procedure here.


1. Note the lens focus ring markings “NEAR “ and “FAR.”




- Place a target with high resolution features in the center of the field of view. Most business cards are adequate for this step. You may leave the door open and focus by ambient light.



- Launch the software, click on the acquire button, then zoom in completely in the preview window. The zoom option will grey out when you reach the end of the range.



- Click on the Auto Expose option.



- Loosen the focus ring locking screw.



- Adjust focus until image is clearly focused then turn focus ring until the image is slightly out of focus in the "NEAR" direction.




- 7. Slowly adjust focus ring only in the “FAR” direction until the image is perfectly focused.




- 8. Lock focus ring in place with focus ring locking screw.



- 9. Remove target.









- 10. Zoom out completely (tool will grey out).
- 11. Take a test image to verify focus.





## Appendix C—Regulatory

### Waste Electrical and Electronic Equipment (WEEE)

ENG		This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.
FRA		Ce symbole indique que les déchets relatifs à l'équipement électrique et électronique ne doivent pas être jetés comme les ordures ménagères non-triées et doivent être collectés séparément. Contactez un représentant agréé du fabricant pour obtenir des informations sur la mise au rebut de votre équipement.
GER		Dieses Symbol kennzeichnet elektrische und elektronische Geräte, die nicht mit dem gewöhnlichen, unsortierten Hausmüll entsorgt werden dürfen, sondern separat behandelt werden müssen. Bitte nehmen Sie Kontakt mit einem autorisierten Beauftragten des Herstellers auf, um Informationen hinsichtlich der Entsorgung Ihres Gerätes zu erhalten.
ITA		Questo simbolo indica che i rifiuti derivanti da apparecchiature elettriche ed elettroniche non devono essere smaltiti come rifiuti municipali indifferenziati e devono invece essere raccolti separatamente. Per informazioni relative alle modalità di smantellamento delle apparecchiature fuori uso, contattare un rappresentante autorizzato del fabbricante.
SPA		Este símbolo indica que el equipo eléctrico y electrónico no debe tirarse con los desechos domésticos y debe tratarse por separado. Contacte con el representante local del fabricante para obtener más información sobre la forma de desechar el equipo.
SWE		Denna symbol anger att elektriska och elektroniska utrustningar inte får avyttras som osorterat hushållsavfall och måste samlas in separat. Var god kontakta en auktoriserad tillverkarrepresentant för information angående avyttring av utrustningen.

### CE Conformity

The following Omega Fluor Imaging Systems,

Models: **Omega Fluor** and **Omega Fluor Plus**

Are in conformity with the provisions of the following EC Directives, including all amendments, and national legislation implementing these directives:

Low Voltage Directive 2006/95/EC

EMC Directive 2004/108/EC

And that the following harmonized standards have been applied:

EN61010-1: 2001

EN61326-1: 1997+A1:1998+A2:2001+A3:2003 Class A

EN61000-4-2, EN61000-4-3, EN61000-4-4, EN61000-4-5, EN61000-4-6, EN61000-4-11

**Protection category: IP20 according IEC 60529**

## Copyright and Trademark Information

All goods and services are sold subject to the terms and conditions of sale of the company within Aplegen International Inc. which supplies them. Aplegen, Inc. reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact an Aplegen, Inc. representative for the current information.

### Aplegen, Inc.

5880 West Las Positas Blvd., Suite 34  
Pleasanton, CA 94588-8552  
USA

Phone: (925) 225-2100  
Fax: (925) 463-3416  
Email: info@aplegen.com  
Tech Support: (925) 463-3410  
Tech Support Fax: (925) 463-3416  
Email: support@aplegen.com

## Safety Warnings

### UV Safety Precautions



The Omega Fluor and Omega Fluor Plus systems come with a built-in Ultra-Violet (UV 200-400nm) Trans-illuminator inside the system. Exposure to UV radiation can cause permanent damage to the eyes and skin. The system enclosure confines the radiation within the system and shields the user from exposure. The system is also equipped with a two-way safety interlock switch which automatically cuts off the power to the trans-illuminator when the door is open during normal use.

The Omega Fluor and Omega Fluor Plus imaging systems belong to Class A equipment, and fulfills the limit values of table 3 but not table 4 of EN 61326:1997+A1+A2+A3.

It may become necessary to defeat the safety lock or operate the transilluminator outside the system for service. In this case, be sure to use the following safety precautions:

- Always wear UV-protected eyewear that is specified by the manufacturer as providing protection at the wavelength(s) used, making sure that the eyewear protects any areas wear radiation may come through (UV sunglasses may not prevent UV radiation from coming in through the sides or around the lenses).
- Always cover all skin that may be exposed to UV light, especially the face, neck, hands, and arms.
- Always make sure that any UV protection devices (such as the safety switch on the light cabinet apparatus) are working properly. If not, discontinue use until the device(s) are properly repaired.
- Please use only UV lamps in the trans-illuminator.

## Electrical Safety Precautions

Be sure to take proper precautions when handling any electrical equipment. NEVER work on any live circuit, fixture, receptacle, or switch. Safety rules you should follow whenever working with any electrical appliance include:

- Always shut off power at the main disconnect before changing a fuse.
- Always shut off power to the circuit before repairing or replacing a switch, receptacle, or fixture.
- Always tape over the main switch, empty fuse socket, or circuit breaker you are working on.
- Always check that the circuit is dead before beginning work on it. Using a circuit tester or voltmeter can help you determine this.
- Always unplug any appliance before repairing it.



The earth terminal, intended for connection to an external protective conductor for protection against electric shock in case of a fault, is located on the inside of the back panel.

### Protective earth terminal

## Hot Surface Warning



Under normal conditions, the temperature of glass surface of UV transilluminator is below 50 °C and safe to touch. However, if the system malfunctions, it is possible that the glass surface temperature exceeds 80°C. Please exercise caution when touch the glass surface with hand when this occurs.

## System Specifications


Camera resolution	5 MP
EPI-illumination	White LED
Trans-Illumination	302 nm/365nm
Lens	8mm/F1.4, ultra low distortion
Maximum Field of View	20 cm x 24 cm
Image output	16-bit Tiff, Jpeg
Power requirement	DC 12V 5A < 60 W
Dimensions (W x D x H)	34.6cm x 31.1cm x 68.6cm
Weight	21.3 kg

## Appendix D—Mitsubishi Electric P-95D Printer Driver Installation


Windows 8 Driver Signature Enforcement is a new feature of Windows 8 that has certain signature requirements not yet implemented in the Mitsubishi Electric P-95D Printer. In order to install the printer driver you must disable the Driver Signature Enforcement on startup, install the driver, and then reboot the system. The Driver Signature Enforcement will reset upon rebooting, but once the driver has been installed it will remain. You will need a USB keyboard to complete the driver installation.

To install the Mitsubishi Electric P-95D Printer Driver please take the following steps:

- Open the Charms Bar, swipe with your finger from the far right side of the screen, to the left.
- Select **Settings**.



- Select **Change PC Settings**.




- Select **General**.




- Scroll down the list, Under Advanced Startup select **Restart Now**. The system will restart.




- Select **Troubleshoot**.




- Select **Advanced Options**.




- Select **Startup Settings**.



- Select **Restart**.




- Use the USB keyboard to select number 7 (hit **F7** or the number 7 key). The system will restart. You are now ready to install the printer driver.



- To start, plug in the printer USB.
- Turn on the printer and follow the directions that appear on the screen.

- Or, you may install the printer from the control panel. To reach the control panel, open the Charms bar, select **Control Panel**.

Detailed driver installation instructions are available in the Mitsubishi Electric P-95D Printer Installation Materials.



## Appendix E—System Requirements

The following are the minimum system requirements to run the Omega Fluor system software. Aplegen recommends using a system with higher than the minimum requirements for optimum performance.


- Windows XP, SP3 and Windows 7, Home or Pro Editions
- .net Framework 3.5 or above
- Pentium 300-megahertz (MHz) processor or faster
- 128 megabytes (MB) of RAM or greater
- 1.5 gigabytes (GB) Hard Drive or greater
- At least one available USB port in addition to those needed for the keyboard, mouse, and optional printer.

## Appendix F—Installing or Changing Emission Filters

The Omega Fluor Plus has a three position sliding filter wheel. Optional filters may be purchased from Aplegen and installed by the user for those applications that require an emission band path that is different from the standard Orange filter (590/50nm).

Installation is quick and requires no tools. Aplegen recommends that you wear disposable glove while handling filters to avoid getting fingerprints on the glass. In the case that prints are left behind they can be cleaned with ethanol and lens paper. Do not use regular paper towels as these may scratch the filter.

First, locate an open filter position. When looking up into the cabinet you will see the Orange filter in Position 1. Position 2 is the middle position, and Position 3 is on the right hand side as you are looking into the cabinet.



Remove the lens from the packaging and insert it into the open filter position. The filter ring is threaded. Turn the filter clockwise to tighten. Take care not to cross the threads.



## Appendix G—USB Hub Installation Guide

The Omega Fluor Plus comes with an optional USB hub for those systems that are not going to be connected to a wireless network. Please use the following guide for easy installation of the hub.

1. Affix the USB hub to a convenient location on the cabinet. The Hub includes double stick tape for placement.



2. Remove the camera cable from the rear cubby.



3. Disconnect the camera USB from the right angle adapter.



4. Slide the camera USB cable through the cable guide and plug it into the Hub.



5. Connect the Hub USB mini cable to the Hub.



6. Connect the Hub USB mini cable to the right angle adapter.



7. Plug in the USB Hub power supply.



8. Roll up excess USB cable and tuck it back into the rear cubby.



## Appendix H—Basic Tablet Operations

The Omega Fluor Plus a touch screen tablet PC to run the imaging station. For in depth directions refer to the Acer Tablet User Manual. However, here are some basic navigation tips to help you get started.

- Use your finger as you would a mouse.
- One touch or tap to select.
- To Scroll tap once, then drag the scroll bar as desired.
- Double tap an icon to launch an application.
- To “right click” touch and hold your finger to the screen. A circle will appear around the point, when it completes, lift your finger and the options menu will appear.
- To exit the menu, tap anywhere on the desktop.
- To enter text information touch and hold on the far left side of the screen. This will cause the virtual keyboard to appear as a tab. Touch the tab to select.
- The windows button on the lower left side of the tablet can be used to open the Windows Main Menu. You can use the shortcut Windows logo key (+D). Please see Omega Fluor Windows 8 Software Installation to see the icon.
- For ease of installation printer drivers for the Mitsubishi P93D Thermal Printer have been pre-installed. To complete the installation simply plug in the USB cable and follow the prompts on the screen.