

## Equilibration of 24 cm IPG Strips for DALT 2D Gel

### Always:

- ◆ Wear powder-free disposable gloves.

### Equilibration

1. Prepare two equilibration solutions in gelcompany IPG-strip equilibration buffer:

**DTT solution:** Weigh urea and dithiothreitol (DTT), and add the equilibration buffer according to table 1 (below) and dissolve completely.

**IAA solution:** Weigh urea and iodoacetamide (IAA), and add the equilibration buffer according to the table 1 and dissolve completely.

2. Using gelcompany equilibrator (fig 1. cat # 1003-04) on an orbital shaker (30 rev/min) equilibrate each strip in DDT solution (6 mL) for 15 min.

Move the strip to a slot containing IAA solution (6 mL) and equilibrate for a further 15 min. Discard the solutions after the 2<sup>nd</sup> equilibration.




Fig. 1

*This protocol forms part of a series that describe the use of gelcompany products in the 2D-gel workflow. For other protocols visit our website.*

| Number of strips | Urea [g] | DTT [mg] | IAA [mg] | E. Buffer [mL] | Total volume [mL] |
|------------------|----------|----------|----------|----------------|-------------------|
| 1                | 1.8      | 50       | -        | 5              | 6                 |
| 1                | 1.8      | -        | 125      | 5              | 6                 |
| 2                | 3.6      | 100      | -        | 10             | 12                |
| 2                | 3.6      | -        | 250      | 10             | 12                |
| 3                | 5.4      | 150      | -        | 15             | 18                |
| 3                | 5.4      | -        | 375      | 15             | 18                |
| 4                | 7.2      | 200      | -        | 20             | 24                |
| 4                | 7.2      | -        | 500      | 20             | 24                |
| 5                | 9        | 250      | -        | 25             | 30                |
| 5                | 9        | -        | 625      | 25             | 30                |
| 6                | 10.8     | 300      | -        | 30             | 36                |
| 6                | 10.8     | -        | 750      | 30             | 36                |
| 7                | 12.6     | 350      | -        | 35             | 42                |
| 7                | 12.6     | -        | 875      | 35             | 42                |
| 8                | 14.4     | 400      | -        | 40             | 48                |
| 8                | 14.4     | -        | 1,000    | 40             | 48                |
| 9                | 16.2     | 450      | -        | 45             | 54                |
| 9                | 16.2     | -        | 1,125    | 45             | 54                |
| 10               | 18       | 500      | -        | 50             | 60                |
| 10               | 18       | -        | 1,250    | 50             | 60                |
| 11               | 19.8     | 550      | -        | 55             | 66                |
| 11               | 19.8     | -        | 1,375    | 55             | 66                |
| 12               | 21.6     | 600      | -        | 60             | 72                |
| 12               | 21.6     | -        | 1,500    | 60             | 72                |

Table 1. Preparing the equilibration buffers