

PROTOCOL

Thermocycle Testing Kit



Part Number: 20K-100
Batch: 1007-1
Storage: Store @ -20°C
Batch Number: Marked on tube

Tube 1: OK Mix

Volume: 1.0 ml
Color of Lid: Yellow

Tube 2: DNA/Primers

Volume: 1.0 ml
Color of Lid: Green

Protocol:

1. Transfer 10 μ l of OK mix (**Tube 1**) into the PCR tube. Add 10 μ l DNA/Primers (**Tube 2**).
2. Overlay with mineral oil if necessary.
3. Place in a Thermal Cycler. If possible, set Thermal Cycler to **ramp at medium rate** (1.2 to 2.5°C/sec)

Cycling Profile:

Initial denaturation step: 95°C for 5 mins


Then cycle 31 times:

Step 1: 95°C for 30 secs

Step 2: 66°C for 60 secs

Step 3: 72°C for 60 secs

After cycling, load 10 μ l onto 1.7% agarose gel and electrophorese alongside a 100 bp DNA ladder.



Expected fragment sizes: 360 bp (bottom), 550 bp (middle) and 650 bp (top).

The 350 bp fragment should be the brightest, then the 550 and then the 650 bp.

G = Good machine: Top band just visible. Machine working correctly.

B = Bad machine: Middle and top band too strong.
Machine cooling below the required annealing temp.

U = Ugly machine: No top band and middle band very faint.
Machine not cooling to the required annealing temp.