



96 ASSAYS ONE PLATE

FOR MULTIPLE ASSAYS



ArrayPlate 96




- 96 (7mm x 7mm) printing surfaces
- round openings
- accommodates septamat

CAT# AHT96-01

assembly includes lower carrier tray, upper structure plate, plate, printing template, septamat and gasket (patent pending)

ArraySlide


- uses standard slides
- easy loading by multi-channel pipettes or automated robots
- utilizes, leak free, compression fit gasket
- no contaminating adhesives
- embossed with letter and number panels for easy identification for each subarray
- durable construction
- will not break slides



ACCESSORIES



ArraySlide 16



ArraySlide 16-4



ArraySlide 24



ArraySlide 24-4



ArrayPlate 96

ArraySlide 16 & ArraySlide 16-4

Gasket	CAT# AUA16
Glass Slide, uncoated	CAT# GEP31-12
Glass Slide, epoxy coated	CAT# AEP16-05

ArraySlide 24 & ArraySlide 24-4

Gasket, for CAT# AHF24	CAT# AUA24
Gasket, for CAT# AHS24-4	CAT# AUA24-4S
Gasket, for CAT# AHM24-4	CAT# AUA24-4M
Glass Slide, uncoated	CAT# GEP31-12
Glass Slide, epoxy coated	CAT# AEP24-05
24 Septa Mat	CAT# SEP24-4
96 Septa Mat	CAT# SEP96-20


ArrayPlate 96

Gasket	CAT# AUA96
Glass Plate, uncoated	CAT# GEP96-1
Glass Plate, epoxy coated	CAT# AEP96-05
96 Septa Mat	CAT# SEP96-20



16 ASSAYS

ONE SLIDE OR FOUR SLIDES




ArraySlide 16

- 16 (7mm x 7mm) printing surfaces
- rectangular openings

CAT# AHF16

assembly includes lower carrier tray, upper structure plate, slide, printing template and gasket (patent pending)




ArraySlide 16-4

- 64 (7mm x 7mm) printing surfaces
- rectangular openings

CAT# AHS16-4


assembly includes lower carrier tray, upper structure plate, slides, printing template and gasket (patent pending)



665 Third Street, Suite 240
San Francisco, CA 94107
415-247-8760 (tel)
415-247-8765 (fax)
www.gelcompany.com

24 ASSAYS

ONE SLIDE OR FOUR SLIDES



ArraySlide 24

- 24 (6mm x 7mm) printing surfaces
- round openings
- accommodates septamat

CAT# AHF24

assembly includes lower carrier tray, upper structure plate, slide, printing template, septamat and gasket (patent pending)



ArraySlide 24-4

- 96 (6mm x 7mm) printing surfaces

CAT# AHS24-4

-rectangular openings
assembly includes lower carrier tray, upper structure plate, slides, printing template and gasket (patent pending)

CAT# AHM24-4

-round openings
-accommodates septamat
assembly includes lower carrier tray, upper structure plate, slide, printing template, septamat and gasket (patent pending)

