

A BETTER WAY TO EXCISE BANDS FROM GELS

GeneCatcher

Disposable gel excision tips

- Designed for removing bands from agarose gels
- One handed operation
- Push button gel release
- Push button tip release
- Avoid cross contamination
- Tips fit standard 1000 μ l pipettors



GeneCatcher Tips

CAT# PKB4.0 - 250/bag

4.0mm x 1mm

CAT# PKB6.5 - 250/bag

6.5mm x 1mm




GeneCatcher Tips - Racks

CAT# PKB4.0-R


5 racks of 48

CAT# PKB6.5-R

5 racks of 48



ACCESSORIES



ISOLATOR

agarose for gel electrophoresis

CAT# AGD050 - 50g

CAT# AGD100 - 100g

CAT# AGD250 - 250g



665 Third Street, Suite 240
San Francisco, CA 94107
415-247-8760 (tel) | 415-247-8765 (fax)
www.gelcompany.com

ARE YOU PAYING TOO MUCH FOR ABI GENETIC ANALYZER ACCESSORIES?



Regenerate your used array for a fraction of the cost of a new one. Compatible with all array types.



Detection plates. Optically clear for real time PCR. Barcoding available. 96 and 384 well formats.



*10X CE Buffer with EDTA
Each lot is function tested.*

ABI310

	ABI#	SIZE	GELCO#
Septa for .5 mL Tubes	401956	500pkg	SEP05-500
Sample Tubes (0.5mL)	401957	500pkg	SET500
10X CE Buffer with EDTA	402824	25mL	DAD-025
310 Capillaries, 47cm	402839	5pkg	IAD47-5
310 Capillaries, 61cm	402840	2pkg	IAD61-2
500 ROX Ladder	401734	800rxn	MRK500
1000 ROX Ladder	n/a	800rxn	MRK1000
BetterBuffer BigDye Enhance	n/a	3mL	DAF-10
Sequencing Grade Water	n/a	25mL	DWW-025
Deionized Formamide	4311320	25mL	DFM-025

ABI3100/3130XL

	ABI#	SIZE	GELCO#
Regenerated Array	4315930	1ea	AR3100
10X CE Buffer with EDTA	402824	25mL	DAX-025
96-Well Plate Septa	4315933	20pkg	SEP96-20
96-Well Detection Plates	4316813	25pkg	MAP100C
500 ROX Ladder	401734	800rxn	MRK500
1000 ROX Ladder	n/a	800rxn	MRK1000
BetterBuffer BigDye Enhance	n/a	3mL	DAF-10
Sequencing Grade Water	n/a	25mL	DWW-025
Deionized Formamide	4311320	25mL	DFM-025

ABI3730/3730XL

	ABI#	SIZE	GELCO#
Regenerated Array	4339252	1ea	AR3730
10X CE Buffer with EDTA	4335613	1L	DAB-01
Heat Seal Film Roll	4337570	1pkg	GC3730
96-Well Plate Septa	4315933	20pkg	SEP96-20
96-Well Detection Plates	4316813	25pkg	MAP100C
500 ROX Ladder	401734	800rxn	MRK500
1000 ROX Ladder	n/a	800rxn	MRK1000
BetterBuffer BigDye Enhance	n/a	3mL	DAF-10
Sequencing Grade Water	n/a	25mL	DWW-025
Deionized Formamide	4311320	25mL	DFM-025
Heat Seal Removing Tool	n/a	1ea	LSP100



665 Third Street, Suite 240
San Francisco, CA 94107
415-247-8760 (tel) | 415-247-8765 (fax)
www.gelcompany.com