

NEON PIPETTORS IN 5 FLAVORS



- accepts universal standard tips
- contoured finger holds
- bright colors
- easy turn volume adjustment from 2 places
- lightweight (under 3oz.)

custom orders welcome...
our quality is guaranteed.

1-800-256-8596

gel
COMPANY

665 Third Street, Suite 240 | San Francisco, CA 94107
415-247-8760 (tel) | 415-247-8765 (fax)
www.gelcompany.com

Green P2: 0.1µl to 2.0µl adjustable volume

precision: $\pm 5\%$ @ 2µl

QDLP02 <i>standard length pipette</i>	\$109
QDSP02 <i>mini pipette</i>	\$109
TP02-R <i>10 racks of 96 tips</i>	\$34
TP02BULK <i>1000 tips</i>	\$14

Pink P10: 1µl to 20µl adjustable volume

precision: $\pm 1\%$ @ 10µl, $\pm 2.5\%$ @ 2µl,

QDLP10 <i>standard length pipette</i>	\$109
QDSP10 <i>mini pipette</i>	\$109
TP10-R <i>10 racks of 96 tips</i>	\$34
TP10BULK <i>1000 tips</i>	\$14

Purple P20: 0.1µl to 20µl adjustable volume

precision: $\pm 1\%$ @ 20µl, $\pm 5\%$ @ 2µl

QDLP20 <i>standard length pipette</i>	\$109
QDSP20 <i>mini pipette</i>	\$109
TP20Y-R <i>10 racks of 96 tips</i>	\$34
TP20YBULK <i>1000 tips</i>	\$14

Orange P200: 10µl to 200µl adjustable volume

precision: $\pm 0.8\%$ @ 200µl, $\pm 2.5\%$ @ 20µl,

QDLP200 <i>standard length pipette</i>	\$109
QDSP200 <i>mini pipette</i>	\$109
TP200Y-R <i>10 racks of 96 tips</i>	\$34
TP200YBULK <i>1000 tips</i>	\$14

Blue P1000: 100µl to 1000µl adjustable volume

precision: $\pm 0.8\%$ @ 1000µl, $\pm 3\%$ @ 200µl,

QDLP1000 <i>standard length pipette</i>	\$109
QDSP1000 <i>mini pipette</i>	\$109
TP1000B-R <i>10 racks of 96 tips</i>	\$34
TP1000BULK <i>1000 tips</i>	\$15

Set of all FIVE pipettors

QDLPSET <i>standard length pipettors</i>	\$439
QDSPSET <i>mini pipettors</i>	\$439