

BRIGHTEN UP YOUR LAB WITH A NEON LAMP MICROFUGE

Crystal Eight

Transparent microfuge with bright NEON lamp

- attractive design
- holds eight tubes (1.5ml to 2.2ml)
- also holds 0.4ml to 0.5ml PCR tubes
- powerful, quiet drive motor
- 10,000 RPM's, 5.6G's
- built in auto stop

Crystal Eight Fuge (blue)

CAT# VS-100B

Crystal Eight Fuge (green)

CAT# VS-100G

Crystal Eight Fuge (red)

CAT# VS-100R



665 Third Street, Suite 240
San Francisco, CA 94107
415-247-8760 (tel)
415-247-8765 (fax)
www.gelcompany.com