



HPE POWER SUPPLY 1500 TECHNICAL SPECIFICATION

Catalogue Number

PPP-001-TI

Control Mode

Constant voltage
 Constant current
 Constant power
 Automatic cross-over
 Automatic programming

Output

Voltage	30 V - 1500 V
Current	0.5 mA - 400 mA
Power	1 W - 300 W

Program steps

Voltage	1 V
Current	1 mA
Power	0.1 W
Timer	1 min/1 h
Volt hour	1 Vh

AC/DC Range

	200 V - 260 V AC	50/60 Hz
Optional	90 V - 130 V AC	50/60 Hz

Power consumption

max. 420 W

Environmental Working range

Temperature	5 C - 40 °C
Humidity	0 % - 95 %
nicht kondensierend	

Dimension

30 × 33 × 10.8 cm

Weight

5.46 Kg