

MATERIAL SAFETY DATA SHEET

Goldmine 10X CE running buffer



Part Number: DAB-XX
Version: 2.0
Storage: Store @ +4.0°C or room temperature for 18 months
Related Accessories: MAP900 96 well microtiter plates
DAF-10 BigDye dilution buffer
2mcl-10 microCLEAN PCR clean up

SAFETY DATA SHEET

DATE PRODUCED: April 2009

- CHEMICAL IDENTIFICATION:** DAB-XX, 10X CE running buffer, Blend
- COMPOSITION/INFORMATION ON INGREDIENTS:** not hazardous
- HAZARDS IDENTIFICATION:**
LABEL PRECAUTIONARY STATEMENTS: irritating to eyes, respiratory system and skin.
- FIRST-AID MEASURES:**
In case of contact with eyes, immediately flush eyes with copious amounts of water for at least 15 mins;
In case of contact with skin, immediately wash skin with soap and copious amounts of water; If inhaled, move to fresh air and monitor breathing; If not breathing give artificial respiration; If breathing is difficult, give oxygen; If swallowed, wash out mouth with water provided person is conscious; Call a physician.
- HANDLING AND STORAGE:**
This product should be handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice.
Storage conditions at +4°C or room temperature.
- PHYSICAL AND CHEMICAL PROPERTIES:** Colorless liquid
- REACTIVITY:**
INCOMPATIBILITIES: Strong oxidising agents
HAZARDOUS COMBUSTION OR DECOMPOSITION PRODUCTS:
Toxic fumes of: Carbon Monoxide, Carbon Dioxide, Hydrogen Chloride.
HAZARDOUS POLYMERISATION: Will not occur
- TOXICOLOGICAL INFORMATION:** Harmful if swallowed
- TRANSPORT INFORMATION:** Non-hazardous for road, sea and air freight
- REGULATORY INFORMATION:** Irritant (R 36/37/38)
Irritating to eyes, respiratory system and skin (S 26)
In case of contact with eyes, rinse immediately with copious amounts of water and seek medical advice (S 36)
Wear suitable protective clothing.
- ECOLOGICAL INFORMATION:** Data not available

MDAB-09