


# OPTICALLY CLEAR 2D CASSETTES 30-40% OFF FOR ETTAN DALT® OR BIO-RAD®

## Low Fluorescent Cassettes for 2-D Gels Compatible with Amersham® DALT and Bio-Rad PROTEAN® Systems



- precision cut glass
- available with or without durable silicone hinge
- bonded spacers (sticky spacers for Ettan DALT; glass spacers for Bio-Rad)
- available in 1.0mm, 1.5mm and 2.0mm spacer thickness

### Catalog #

GHY2620-100

GHY2620-150

EDF06

EDF16

GBD2020-100

GBD2020-150

GBD2520-100

GBD2520-150

GBN2020-100

GBN2020-150

GBN2520-100

GBN2520-150

RLB-10

### Description

Optically clear Ettan DALT cassette, no hinge, 1.0mm

Optically clear Ettan DALT cassette, no hinge, 1.5mm

Ettan DALT filler sheets, 1.0mm, 6 units

Ettan DALT separator sheets, 0.5mm, 16 units

Bio-Rad Dodeca PLUS cassette, hinged, 20 x 20.5cm, 1.0mm

Bio-Rad Dodeca PLUS cassette, hinged, 20 x 20.5cm, 1.5mm

Bio-Rad Dodeca PLUS cassette, hinged, 25 x 20.5cm, 1.0mm

Bio-Rad Dodeca PLUS cassette, hinged, 25 x 20.5cm, 1.5mm

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 20 x 20.5cm, 1.0mm

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 20 x 20.5cm, 1.5mm

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 25 x 20.5cm, 1.0mm

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 25 x 20.5cm, 1.5mm

Safeplate Rack, 27 x 15 x 19cm, black ABS, 10 slot

Many more sizes and options available. Please see our website for more information.  
Accessories include 2D gel spotpickers, lane cutters and filler sheets for gel casting units.



665 Third Street, Suite 240

San Francisco, CA 94107

415-247-8760 (tel)

415-247-8765 (fax)

[www.gelcompany.com](http://www.gelcompany.com)

*Amersham and ETTAN are trademarks of GE Healthcare.*

*Bio-Rad and PROTEAN are registered trademarks of Bio-Rad Laboratories.*

**VISIT OUR WEBSITE FOR NEW PRODUCTS AND LAB SPECIALS**

**“your discount accessory supplier”**

# A BETTER WAY TO EXCISE SPOTS & BANDS FROM GELS

## OneTouch

2-D gel spot picker



- no more slicing with scalpels or blades
- for removal from polyacrylamide and 6mm agarose gels
- one motion removes spot
- push button gel release
- ergonomic and time saving
- available in 1.5mm, 3.0 and 5.0 diameters

**OneTouch - 1.5mm diameter - 1 unit**

**CAT# P2D1.5**

**OneTouch - 3.0mm diameter - 1 unit**

**CAT# P2D3.0**

**OneTouch - 5.0mm diameter - 1 unit**

**CAT# P2D5.0**



## G-Tubes

Low protein binding microtubes for in-gel digestion

**G-Tube - 0.65ml - 100/pack**

**CAT# LPB65-100**

**G-Tube - 0.65ml - 50/pack**

*individually wrapped*

**CAT# LPBT65-50**

## PowerPicker

1D gel band remover




- thin ramped edge for clean, efficient gel band removal
- 5.0mm x 1.0mm rectangular cutter
- surgical grade stainless steel
- comfortable push grip

**Rectangular Band Picker - 1 unit**

**CAT# PDB5.0**

## OneTouch GridCutter

Multiple band excision



- will excise and remove 40 bands
- each 4mm x 1mm
- disposable blades and mount are sold separately
- mount comes with special angled forceps for band removal

**OneTouch GridCutter - 10/bag**

**CAT# MEE-1x5**

**OneTouch Mount - 1 unit**

**CAT# MEF-1.5**



665 Third Street, Suite 240 | San Francisco, CA 94107

415-247-8760 (tel) | 415-247-8765 (fax)

[www.gelcompany.com](http://www.gelcompany.com)

BGHY2620-100-09