


OPTICALLY CLEAR 2D CASSETTES 30-40% OFF FOR ETTAN DALT® OR BIO-RAD®

Low Fluorescent Cassettes for 2-D Gels Compatible with Amersham® DALT and Bio-Rad PROTEAN® Systems



- Precision cut glass
- Available with or without durable silicone hinge
- Bonded spacers (sticky spacers for Ettan DALT; glass spacers for Bio-Rad)
- Available in 1.0mm, 1.5mm and 2.0mm spacer thickness

Catalog

Description

GHY2620-100

Optically clear Ettan DALT cassette, no hinge, 1.0mm

GHY2620-150

Optically clear Ettan DALT cassette, no hinge, 1.5mm

EDF06

Ettan DALT filler sheets, 1.0mm, 6 units

EDF16

Ettan DALT separator sheets, 0.5mm, 16 units

GBD2020-100

Bio-Rad Dodeca PLUS cassette, hinged, 20 x 20.5cm, 1.0mm

GBD2020-150

Bio-Rad Dodeca PLUS cassette, hinged, 20 x 20.5cm, 1.5mm

GBD2520-100

Bio-Rad Dodeca PLUS cassette, hinged, 25 x 20.5cm, 1.0mm

GBD2520-150

Bio-Rad Dodeca PLUS cassette, hinged, 25 x 20.5cm, 1.5mm

GBN2020-100

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 20 x 20.5cm, 1.0mm

GBN2020-150

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 20 x 20.5cm, 1.5mm

GBN2520-100

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 25 x 20.5cm, 1.0mm

GBN2520-150

Optically clear Bio-Rad Dodeca PLUS cassette, no hinge, 25 x 20.5cm, 1.5mm

RLB-10

Safeplate Rack, 27 x 15 x 19cm, black ABS, 10 slot

Many more sizes and options available. Please see our website for more information.
Accessories include 2D gel spotpickers, lane cutters and filler sheets for gel casting units.

Amersham and ETTAN are trademarks of GE Healthcare.

Bio-Rad and PROTEAN are registered trademarks of Bio-Rad Laboratories.



665 Third Street, Suite 240
San Francisco, CA 94107
415-247-8760 (tel) | 415-247-8765 (fax)
www.gelcompany.com

VISIT OUR WEBSITE FOR NEW PRODUCTS AND LAB SPECIALS

“your discount accessory supplier”

A BETTER WAY TO EXCISE SPOTS & BANDS FROM GELS

OneTouch

2-D gel spot picker



- No more slicing with scalpels or blades
- For removal from polyacrylamide and 6mm agarose gels
- One motion removes spot
- Push button gel release
- Ergonomic and time saving
- Available in 1.5mm and 3.0mm diameters

OneTouch - 1.5mm diameter - 1 unit

CAT# P2D1.5

OneTouch - 3.0mm diameter - 1 unit

CAT# P2D3.0

PowerPicker

1D gel band remover




- Thin ramped edge for clean, efficient gel band removal
- 5.0mm x 1.0mm rectangular cutter
- Surgical grade stainless steel
- Comfortable push grip

Rectangular Band Picker - 1 unit

CAT# PDB5.0

OneTouch GridCutter

Multiple band excision



- Will excise and remove 40 bands
- Each 4mm x 1mm
- Disposable blades and mount are sold separately
- Mount comes with special pick for band removal

OneTouch GridCutter - 10/bag

CAT# MEE-1x5

OneTouch Mount - 1 unit

CAT# MEF-1.5



G-Tubes

Low protein binding microtubes for in-gel digestion

G-Tube - 0.65ml - 100/pack

CAT# LPB65-100

G-Tube - 0.65ml - 50/pack

individually wrapped

CAT# LPBT65-50



665 Third Street, Suite 240
San Francisco, CA 94107
415-247-8760 (tel) | 415-247-8765 (fax)
www.gelcompany.com