

ARE YOU PAYING TOO MUCH FOR DEEP WELL PLATES?

ARE YOU PAYING TOO MUCH FOR PCR PLATES?



96 well deep well plate

- Thin walled
- 1.2ml max well volume (when sealed)
- Cut-away corner for precise orientation
- Compatible with ABI's SOLiD™ system

96w Deep Well Plate

CAT# MAP471C
50 units - Clear





OPTICALLY CLEAR SEALING FILM AVAILABLE

Optically transparent polyester film with a strong, ultra-smooth, non-absorbing, non-flourescing medical-grade adhesive for superior performance in real-time PCR application. Brilliant ultra-high optical clarity. Available in 2 mil. thickness and for more rigid, 5 mil. thickness. Certified DNase-, RNase-, and nucleaic-acid-free.


CAT# MAF800	ThermalSeal RT, 2 mil thickness	100 units
CAT# MAF850	ThermalSeal RT, 5 mil thickness	100 units

NEW PCR Plates

We offer a variety of NEW PCR plates that are compatible with Applied Biosystems, Biometra, Bio-Rad, Eppendorf, MJ Research and Stratagene cyclers.



FOR TIMES WHEN YOU NEED LESS THAN 96 SAMPLES



96 well non-skirted break away plate

- Can be sectioned into 8 well strips
- 0.2ml max well volume (when sealed)
- Cut-away corner for precise orientation
- Compatible with most major PCR and QPCR Cyclers


96w Non-Skirted, Break Away Plate

CAT# MAP980C
25 units - Clear



665 Third Street, Suite 240
San Francisco, CA 94107
415-247-8760 (tel) | 415-247-8765 (fax)
www.gelcompany.com

LOOKING FOR A BETTER DEAL ON FAST PLATES?



96 well semi-skirted, low profile

- Thin walled
- 0.1ml max well volume
(when used with adhesive and heat seals)
- Cut-away corner for precise orientation
- Clear well for sample visibility
- Designed to be compatible with all ABI Fast Block machines

96w Semi-Skirted, Low Profile Plate

CAT# MAP190C
25 units - Clear