

SIMULTANEOUSLY EXCISE  
MULTIPLE BANDS FROM GELS

A BETTER WAY TO EXCISE  
SPOTS & BANDS FROM GELS


## OneTouch GridCutter

multiple band excision



- Simultaneously excise multiple gel slices
- Disposable blades are numbered for easy identification
- Mount and disposable blades are sold separately
- Mount comes with special pick for band removal
- Available in multiple sizes
- Ask about custom configurations

### ACCESSORIES



#### G-Tubes

Low protein binding microtubes for in-gel digestion

#### CAT# LPB65-100

100/pack - 0.65ml

#### CAT# LPBT65-50

50/pack - 0.65ml  
individually wrapped

#### OneTouch GridCutter

##### CAT# MEE-1X10

Excises 26 bands (1.0mm x 10.0mm)

#### OneTouch Mount

##### CAT# MEF-11

#### OneTouch GridCutter

##### CAT# MEE-1.3X1.3

Excises 36 bands (1.3mm x 1.3mm)

##### CAT# MEE-1X1

Excises 64 bands (1.0mm x 1.0mm)

##### CAT# MEE-2X2

Excises 16 bands (2.0mm x 2.0mm)

#### OneTouch Mount

##### CAT# MEF-8.8

#### OneTouch GridCutter

##### CAT# MEE-1X3.5

Excises 26 bands (1.0mm x 3.5mm)

#### OneTouch Mount

##### CAT# MEF-3.5

#### OneTouch GridCutter

##### CAT# MEE-1X5

Excises 40 bands (1.0mm x 5.0mm)

##### CAT# MEE-2X7

Excises 25 bands (2.0mm x 7.0mm)

#### OneTouch Mount

##### CAT# MEF-1.5

#### OneTouch GridCutter

##### CAT# MEE-1.5X5



Excises 48 bands (1.5mm x 5.0mm)

#### OneTouch Mount

##### CAT# MEF-72

## Proteomics Tools

We offer a variety of proteomics accessories, including OneTouch, PowerPicker, One Touch Plus, ProPicker & OneTouch GridCutter.



**gel**  
company

## REUSABLE SPOT & BAND PICKERS

### OneTouch 2-D gel spot picker



- No more slicing with scalpels or blades
- For removal from polyacrylamide and 6mm agarose gels
- One motion removes spot
- Push button gel release
- Ergonomic and time saving
- Available in 1.5mm and 3.0mm diameters

**CAT# P2D1.5**  
1.5mm diameter

**CAT# P2D3.0**  
3.0mm diameter

### PowerPicker 1-D gel band remover



- Thin ramped edge for clean, efficient gel band removal
- 5.0mm x 1.0mm rectangular cutter
- Surgical grade stainless steel
- Comfortable push grip


**CAT# PDB5.0**  
5.0mm x 1.5mm



665 Third Street, Suite 240 | San Francisco, CA 94107  
415-247-8760 (tel) | 415-247-8765 (fax)  
www.gelcompany.com

## DISPOSABLE SPOT & BAND PICKERS

### OneTouch Plus round gel excision tips



US PATENT 6,702,990

- For removal from 0.75mm, 1.0mm, 1.5mm polyacrylamide and 6mm agarose gels
- One motion removes spot
- No more slicing with scalpels or blades


**OneTouch Plus**  
**CAT# PDM1.5**  
1.5mm diameter

**CAT# PDM3.0**  
3.0mm diameter

**OneTouch Plus Tips**  
**CAT# TSP1.5**  
1.5mm diameter, 10 racks of 96

**CAT# TSP3.0**  
3.0mm diameter, 10 racks of 96

### ProCatcher gel excision tips



- Designed for removing bands from agarose gels
- One handed operation
- Push button gel and tip release
- Avoid cross contamination
- Tips fit standard 1000 µl pipettors

**CAT# PR04.0**  
4.0mm x 1.0mm, 250 tips

**CAT# PR06.5**  
6.5mm x 1.0mm, 250 tips

**CAT# PR04.0-R**  
4.0mm x 1.0mm, 5 racks of 48 tips

**CAT# PR06.5-R**  
6.5mm x 1.0mm, 5 racks of 48 tips