

50bp DNA Ladder RTU (Ready-to-Use)

Cat. No. DDL-012

Size: 50µg / 500µl

Description

A unique combination of PCR products and a number of proprietary plasmids digested with appropriate restriction enzymes to yield 17 fragments, suitable for use as molecular weight standards for agarose gel electrophoresis. The DNA includes fragments ranging from 50-1,500 base pairs. The 200 and 500 base pair bands have increased intensity to serve as reference points. The approximate mass of DNA in each band is provided (0.5 µg a load) for approximating the mass of DNA in comparably intense samples of similar size.

Source

PCR products and double-stranded DNA digested with appropriate restriction enzymes, are phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 1mM EDTA.

Range: 50-1,500 bp

Number of bands: 17

Concentration: 100 µg/ml

Package: 50ug/500ul


Recommended Load: 5 µl / well
Containing orange G as tracking dyes.

Storage

Store at -20°C for 24 months

DNA Mass
(ng/5µl)

Base Pairs



2% TAE agarose gel



www.gelcompany.com