


# MATERIAL SAFETY DATA SHEET

BetterBuffer



Part Number: DAF-XX  
Version: 3.1  
Storage: Store @ -20°C until thawed, then +4.0°C for 4 weeks  
Related Accessories: MAP100 96 well Amp plates  
10X CE running buffers

SAFETY DATA SHEET

DATE PRODUCED: April 2009

1. CHEMICAL IDENTIFICATION: BetterBuffer
2. COMPOSITION/INFORMATION ON INGREDIENTS: Proprietary Mix
3. HAZARDS IDENTIFICATION:  
LABEL PRECAUTIONARY STATEMENTS: irritating to eyes, respiratory system and skin.
4. FIRST-AID MEASURES:  
In case of contact with eyes, immediately flush eyes with copious amounts of water for at least 15 mins.;  
In case of contact with skin, immediately wash skin with soap and copious amounts of water; If inhaled,  
move to fresh air and monitor breathing; If not breathing give artificial respiration; If breathing is difficult,  
give oxygen; If swallowed, wash out mouth with water provided person is conscious; Call a physician.
5. HANDLING AND STORAGE:  
This product should be handled only by those persons who have been trained in laboratory techniques  
and that it is used in accordance with the principles of good laboratory practice.  
Storage conditions at -20°C.
6. PHYSICAL AND CHEMICAL PROPERTIES: Colourless slightly viscous liquid
7. REACTIVITY:  
INCOMPATIBILITIES: Strong oxidising agents  
HAZARDOUS COMBUSTION OR DECOMPOSITION PRODUCTS:  
Toxic fumes of: Carbon Monoxide, Carbon Dioxide, Hydrogen Chloride.  
HAZARDOUS POLYMERISATION: Will not occur
8. TOXICOLOGICAL INFORMATION: Harmful if swallowed
9. TRANSPORT INFORMATION: Non-hazardous for road, sea and air freight
10. REGULATORY INFORMATION:  
Irritant (R 36/37/38)  
Irritating to eyes, respiratory system and skin (S 26)  
In case of contact with eyes, rinse immediately with copious amounts of water and seek medical advice (S 36)  
Wear suitable protective clothing

MDAF-09